





Quantification

Hydrogen peaks on autoclaved samples:
 a) Zry + ZrO₂ b) ZrO₂ (Zry removed chemically)





oxide layer spindled off in hot cell.

Hydrogen peaks on irradiated samples:
 a) Zry + ZrO₂ b) Zry (ZrO₂ removed mechanically)




Neutron Radiography

ring sample
 w. exact geometry.



Voltage Intervals for variable calibration


von einer bestimmten dop



calibration factor k

variable factor v

Calculation of Oxide Thickness



Oxide Thickness s (EC-measurements)

Former Wall Thickness w

Pilling-Bedworth-ratio

Characteristic Data of Samples

Table 4: Characteristic PIE Data of the Investigated Fuel Rods and Samples

Rod/Sample Designation	Liner	Reactor Cycles/ Full Power Days	Burnup [MWd/kg U]	Oxide Thickness [μm]	Diameter [mm]	Wall Thickness [mm]
FR1/H	Zr-alloy 1	5/1656	66.1	66	10.757	0.669
FR2/I	Zr-alloy 1	4/1325	59.4	42	10.747	0.663
FR2/D	Zr-alloy 1	4/1325	59.4	22	10.730	0.676
FR3/G	Zr-alloy 1	1/338	18.9	7	10.672	0.740
FR4/I	Zr-alloy 2	5/1656	66.8	64	10.750	0.675



Results of Fractionated Hydrogen Determination

Results of Fractionated Hydrogen Determination of Fuel Cladding Samples

Sample Designation	Weight [g]	Mass Ratio Oxide/Metal	Fraction	Calculated Weight [g]	Hydrogen [wt.-ppm]
FR1/H-H1/H2	0.4993/0.4373	0.0867	oxide metal	0.0398/0.0348 0.4595	1140/1144 682/702
FR1/H-H3*	0.2758		metal*		651
FR2/I-H1/H2	0.4342/0.4184	0.0554	oxide metal	0.0228/0.0220 0.4114/0.3964	907/1158 415/514
FR2/I-H3*	0.2585		metal*		526
FR2/D-H1/H2	0.4519/0.4089	0.0284	oxide metal	0.0125/0.0113 0.4394/0.3976	821/1125 176/207
FR2/D-H3*	0.3124		metal*		192
FR3/G-H1	0.2999	0.0084	oxide metal	0.0025 0.2974	993 60
FR4/I-H1/H2	0.4208/0.3965	0.0832	oxide metal	0.0323/0.0305 0.3885/0.3660	1077/1008 303/270
FR4/I-H3*	0.3878		metal*		262

*oxide removed by grinding

Results of Fractionated Analysis



DEVELOPING A METHOD FOR FRACTIONATED HYDROGEN DETERMINATION IN FUEL CLADDING SAMPLES

Holger Wiese, Armin Hermann, Raimund Bühner, Michael Steinemann and Gerhard Bart

Laboratory for Materials Behaviour
Paul Scherrer Institute, 5232 Villigen, Switzerland

Introduction

Method

Validating

Application

Conclusions


Introduction

- Hydrogen pick up, an important aspect of cladding corrosion and life limiting criterion of fuel rods for high burnup
- Some hydrogen picked up by ZrO_2
- [H]-determination normally on cladding samples including ZrO_2
- Attempts for separation of [H]-metal vs. [H]-oxide by mechanical removal of ZrO_2
- Development of fractionated hot gas extraction technique at PSI for distinction between [H]-metal and [H]-oxide


Principle of Hydrogen hot gas extraction technique at PSI

- After NDT in hot cells, fuel cladding cut into ring samples, fuel removed, samples UT-cleaned
- Sample transfer to a shielded box with a remote handled LECO-RH-402 Hydrogen Determinator installed
- Cladding rings fed into graphite crucibles
- Crucibles heated to $2100^\circ C$ by RF-induction heating
- Cladding molten together with fluxing agent (Sn)
- N_2 carrier gas used for transportation of hydrogen through heat conductance measuring cell outside of the shielding


Schematic of Modified Hydrogen Analyzer




RF-Power versus (steady state) Crucible Temperature



RF-Power Sweep indicating Hydrogen Pulses from several Reservoirs



Time-Power/Voltage Diagrams P1 and P2



Details of Time-Power Programs P1 and P2

Time [timescale, s]	Step	Power P [%]	
		P1	P2
t_A [-120]	outgassing	100	100
t_r [-30]	stand-by (cooling)	0	0
t_0 [0]*	start, step 1	15	15
t_1 [45]	stop, step 1	15	15
t_2 [46]	start, step 2	90	15
t_3 [136]	end, step 2	90	**
t_2' [166]	switch to P=const.		95
t_3 [186]	end, step 2		95

*start of analysis with sample; **linear power increase up to 95%

Conclusions

- A hot gas extraction method was developed to quantify in one step the fraction of hydrogen present in ZrO_2 layers versus the fraction dissolved/precipitated in the Zry metal phase
- The correlation of the respective hydrogen reservoirs (H_{ZrO_2} and H_{Zry}) to the extraction signals was reached by separate analysis of samples with metal phase only and samples of separated ZrO_2
- Exact measuring of the different masses of hydrogen within the oxide and the metal was achieved by a variable calibration procedure using several certified hydrogen standards
- Water traces as contamination is excluded in the analysis by careful drying of the samples prior to the hot gas extraction process. Moreover the LECO-analyzer is equipped with a water trap of $Mg(ClO_4)_2$

WI43 / 00_frac_H.doc

Conclusions, continued

- The analyzed samples, increasing in burn up and oxide thickness, show a nearly constant (saturated?) hydrogen concentration within the ZrO_2 layer of about 1000 wtppm
- Samples of the same cladding alloy show a constant ratio of mass of hydrogen stored in the oxide to the total mass of hydrogen
- Applying the common hot gas extraction technique without differentiation of H_{ZrO_2} from H_{Zry} leads to a hydrogen concentration overestimation in the metal phase of 5-20%

WI43 / 00_frac_H.doc

Reference:

A. Hermann, H. Wiese, R. Bühner, M. Steinemann, G. Bart; Hydrogen distribution between fuel cladding metal and overlying corrosion layers; Proc. Int. Topical Meeting on Light-Water-Reactor-Performance; April 10-13, 2000; Park City, Utah; 372

WI43 / 00_frac_H.doc

